

Product Datasheet

CD117 antibody (Biotin) (orb1185342)

Description	Rat monoclonal antibody against CD117 conjugated to Biotin
Species/Host	Rat
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Concentration	0.5 mg/mL
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Rat IgG2b, κ
Uniprot ID	P05532
Entrez	16590
Expiration Date	12 months from date of receipt.