

Product Datasheet

CD54 antibody (PerCP/Cy5.5) (orb1185403)

Description	Rat monoclonal antibody against CD54 conjugated to PerCP/Cy5.5
Species/Host	Rat
Reactivity	Mouse
Conjugation	PerCP/Cy5.5
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Concentration	0.2 mg/mL
Storage	Store at 2~8°C. Do not freeze.
Note	For research use only
Isotype	Rat IgG2b, κ
Uniprot ID	P13597
Entrez	15894
Expiration Date	12 months from date of receipt.