

Product Datasheet

CD29 antibody (PerCP/Cy5.5) (orb1185567)

Description	Mouse monoclonal antibody against CD29 conjugated to PerCP/Cy5.5
Species/Host	Mouse
Reactivity	Human
Conjugation	PerCP/Cy5.5
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Mouse IgG1, κ
Uniprot ID	P05556
Entrez	3688
Expiration Date	12 months from date of receipt.