

Product Datasheet

CD45.2 antibody (PE/Cy5) (orb1185680)

Description	Mouse monoclonal antibody against CD45.2 conjugated to PE/Cy5
Species/Host	Mouse
Reactivity	Mouse
Conjugation	PE/Cy5
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Concentration	0.2 mg/mL
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Mouse IgG2a, κ
Uniprot ID	P06800
Entrez	19264
Expiration Date	12 months from date of receipt.