

Product Datasheet

CD11c antibody (PE/Cy7) (orb1185798)

Description	Hamster monoclonal antibody against CD11c conjugated to PE/Cy7
Species/Host	Hamster
Reactivity	Mouse
Conjugation	PE/Cy7
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Concentration	0.2 mg/mL
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Armenian Hamster IgG
Uniprot ID	Q9QXH4
Entrez	16411
Expiration Date	12 months from date of receipt.