

Product Datasheet

CD20 antibody (PE/Cy5.5) (orb1185819)

Description	Mouse monoclonal antibody against CD20 conjugated to PE/Cy5.5
Species/Host	Mouse
Reactivity	Human
Conjugation	PE/Cy5.5
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Mouse IgG2b, κ
Uniprot ID	P11836
Entrez	931
Expiration Date	12 months from date of receipt.