

Product Datasheet

CD4 antibody (orb1185966)

Description	Rat monoclonal antibody to CD4
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	Blocking, Depletion, FC, IHC, IP, WB
Preservatives	0.2 µm filtered in PBS, pH7.2. Azide Free (AF)/Low Endotoxin (LE): Contains no stabilizers or preservatives. Endotoxin level is 2 EU/mg (0.002EU/µg) as Determined by LAL gel clotting assay.
Concentration	0.2 mg/mL
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Rat IgG2b, κ
Uniprot ID	P06332
Entrez	12504
Expiration Date	12 months from date of receipt.