

## Product Datasheet

### CD71 antibody (PerCP/Cy5.5) (orb1185971)

|                            |                                                                               |
|----------------------------|-------------------------------------------------------------------------------|
| <b>Description</b>         | Rat monoclonal antibody against CD71 conjugated to PerCP/Cy5.5                |
| <b>Species/Host</b>        | Rat                                                                           |
| <b>Reactivity</b>          | Mouse                                                                         |
| <b>Conjugation</b>         | PerCP/Cy5.5                                                                   |
| <b>Tested Applications</b> | FC                                                                            |
| <b>Preservatives</b>       | PBS with 0.05% Proclin300, 1% rAlbumin                                        |
| <b>Storage</b>             | Store at 2~8°C and protected from prolonged exposure to light. Do not freeze. |
| <b>Note</b>                | For research use only                                                         |
| <b>Isotype</b>             | Rat IgG2a, $\kappa$                                                           |
| <b>Uniprot ID</b>          | <a href="#">Q62351</a>                                                        |
| <b>Entrez</b>              | <a href="#">22042</a>                                                         |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                               |