

## Product Datasheet

### CD71 antibody (PE/Cy5.5) (orb1185972)

|                            |                                                                               |
|----------------------------|-------------------------------------------------------------------------------|
| <b>Description</b>         | Rat monoclonal antibody against CD71 conjugated to PE/Cy5.5                   |
| <b>Species/Host</b>        | Rat                                                                           |
| <b>Reactivity</b>          | Mouse                                                                         |
| <b>Conjugation</b>         | PE/Cy5.5                                                                      |
| <b>Tested Applications</b> | FC                                                                            |
| <b>Preservatives</b>       | PBS with 0.05% Proclin300, 1% rAlbumin                                        |
| <b>Storage</b>             | Store at 2~8°C and protected from prolonged exposure to light. Do not freeze. |
| <b>Note</b>                | For research use only                                                         |
| <b>Isotype</b>             | Rat IgG2a, κ                                                                  |
| <b>Uniprot ID</b>          | <b>Q62351</b>                                                                 |
| <b>Entrez</b>              | <b>22042</b>                                                                  |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                               |