

Product Datasheet

CD103 antibody (PE) (orb1185977)

Description	Rat monoclonal antibody against CD103 conjugated to PE
Species/Host	Rat
Reactivity	Mouse
Conjugation	PE
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Rat IgG2a, κ
Uniprot ID	Q60677
Entrez	16407
Expiration Date	12 months from date of receipt.