

Product Datasheet

CD161/NK1.1 antibody (PE/Cy7) (orb1186092)

Description	Mouse monoclonal antibody against CD161/NK1.1 conjugated to PE/Cy7
Species/Host	Mouse
Reactivity	Mouse
Conjugation	PE/Cy7
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Concentration	0.2 mg/mL
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Mouse IgG2a, κ
Uniprot ID	P27814
Entrez	17059
Expiration Date	12 months from date of receipt.