

Product Datasheet

CD21 antibody (Biotin) (orb1186149)

Description	Mouse monoclonal antibody against CD21 conjugated to Biotin
Species/Host	Mouse
Reactivity	Human
Conjugation	Biotin
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Concentration	0.5 mg/mL
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Mouse IgG1, κ
Uniprot ID	P20023
Entrez	1380
Expiration Date	12 months from date of receipt.