

Product Datasheet

CD29 antibody (PerCP) (orb1186171)

Description	Mouse monoclonal antibody against CD29 conjugated to PerCP
Species/Host	Mouse
Reactivity	Human
Conjugation	PerCP
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Mouse IgG1, κ
Uniprot ID	P05556
Entrez	3688
Expiration Date	12 months from date of receipt.