

Product Datasheet

IL-4 antibody (APC) (orb1186195)

Description	Rat monoclonal antibody against IL-4 conjugated to APC
Species/Host	Rat
Reactivity	Mouse
Conjugation	APC
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Concentration	0.2 mg/mL
Storage	Store at 2~8°C. Do not freeze.
Note	For research use only
Isotype	Rat IgG1, κ
Uniprot ID	P07750
Entrez	16189
Expiration Date	12 months from date of receipt.