

Product Datasheet

CD120b antibody (PE) (orb1186266)

Description	Hamster monoclonal antibody against CD120b conjugated to PE
Species/Host	Hamster
Reactivity	Mouse
Conjugation	PE
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Concentration	0.2 mg/mL
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Armenian Hamster IgG
Uniprot ID	P25119
Entrez	21938
Expiration Date	12 months from date of receipt.