

Product Datasheet

CD134 antibody (FITC) (orb1186357)

Description	Mouse monoclonal antibody against CD134 conjugated to FITC
Species/Host	Mouse
Reactivity	Human
Conjugation	FITC
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Mouse IgG1, κ
Uniprot ID	P43489
Entrez	7293
Expiration Date	12 months from date of receipt.