

Product Datasheet

HLA-A,B,C antibody (FITC) (orb1186369)

Description	Mouse monoclonal antibody against HLA-A,B,C conjugated to FITC
Species/Host	Mouse
Reactivity	Human
Conjugation	FITC
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Mouse IgG2a, κ
Uniprot ID	P01889
Entrez	3105
Expiration Date	12 months from date of receipt.