

Product Datasheet

IL-4 antibody (orb1186659)

Description	Rat monoclonal antibody to IL-4
Species/Host	Rat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC
Preservatives	PBS with 0.05% Proclin300
Concentration	0.2 mg/mL
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Rat IgG1, κ
Uniprot ID	P05112
Entrez	3565
Expiration Date	12 months from date of receipt.