

Product Datasheet

CD279/PD-1 antibody (FITC) (orb1186726)

Description	Rat monoclonal antibody against CD279/PD-1 conjugated to FITC
Species/Host	Rat
Reactivity	Mouse
Conjugation	FITC
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Concentration	0.5 mg/mL
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Rat IgG2a, κ
Uniprot ID	Q02242
Entrez	18566
Expiration Date	12 months from date of receipt.