

Product Datasheet

CD45RA antibody (PerCP/Cy5.5) (orb1186797)

Description	Mouse monoclonal antibody against CD45RA conjugated to PerCP/Cy5.5
Species/Host	Mouse
Reactivity	Human
Conjugation	PerCP/Cy5.5
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Mouse IgG2b, κ
Uniprot ID	P08575
Entrez	5788
Expiration Date	12 months from date of receipt.