

Product Datasheet

CD62L antibody (PerCP) (orb1186808)

Description	Mouse monoclonal antibody against CD62L conjugated to PerCP
Species/Host	Mouse
Reactivity	Human
Conjugation	PerCP
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Mouse IgG1, κ
Uniprot ID	P14151
Entrez	6402
Expiration Date	12 months from date of receipt.