

Product Datasheet

CD29 antibody (orb1186809)

Description	Mouse monoclonal antibody to CD29
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, IP
Preservatives	PBS with 0.05% Proclin300
Concentration	0.2 mg/mL
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Mouse IgG1, κ
Uniprot ID	P05556
Entrez	3688
Expiration Date	12 months from date of receipt.