

Product Datasheet

CD11c antibody (Biotin) (orb1186848)

Description	Hamster monoclonal antibody against CD11c conjugated to Biotin
Species/Host	Hamster
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Concentration	0.5 mg/mL
Storage	Store at 2~8°C. Do not freeze.
Note	For research use only
Isotype	Armenian Hamster IgG
Uniprot ID	Q9QXH4
Entrez	16411
Expiration Date	12 months from date of receipt.