

Product Datasheet

CD18 antibody (orb1186910)

Description	Rat monoclonal antibody to CD18
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC, IP, WB
Preservatives	PBS with 0.05% Proclin300
Concentration	0.2 mg/mL
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Rat IgG2a, κ
Uniprot ID	P11835
Entrez	16414
Expiration Date	12 months from date of receipt.