

Product Datasheet

CD279/PD-1 antibody (orb1186999)

Description	Mouse monoclonal antibody to CD279/PD-1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC
Preservatives	PBS with 0.05% Proclin300
Concentration	0.2 mg/mL
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Mouse IgG1, κ
Uniprot ID	P20701
Entrez	3683
Expiration Date	12 months from date of receipt.