

Product Datasheet

CD45RO antibody (APC) (orb1187034)

Description	Mouse monoclonal antibody against CD45RO conjugated to APC
Species/Host	Mouse
Reactivity	Human
Conjugation	APC
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Storage	Store at 2~8°C. Do not freeze.
Note	For research use only
Isotype	Mouse IgG2a, κ
Uniprot ID	P08575
Entrez	5788
Expiration Date	12 months from date of receipt.