

Product Datasheet

CD45.2 antibody (PerCP) (orb1187270)

Description	Mouse monoclonal antibody against CD45.2 conjugated to PerCP
Species/Host	Mouse
Reactivity	Mouse
Conjugation	PerCP
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Mouse IgG2a, κ
Uniprot ID	P06800
Entrez	19264
Expiration Date	12 months from date of receipt.