

Product Datasheet

CD54 antibody (PerCP) (orb1187411)

Description	Rat monoclonal antibody against CD54 conjugated to PerCP
Species/Host	Rat
Reactivity	Mouse
Conjugation	PerCP
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Rat IgG2b, κ
Uniprot ID	P13597
Entrez	15894
Expiration Date	12 months from date of receipt.