

## Product Datasheet

### Rat IgM, $\kappa$ Isotype Control (FITC) (orb1187568)

|                            |                                                                               |
|----------------------------|-------------------------------------------------------------------------------|
| <b>Description</b>         | Rat IgM, $\kappa$ Isotype Control (FITC)                                      |
| <b>Species/Host</b>        | Rat                                                                           |
| <b>Conjugation</b>         | FITC                                                                          |
| <b>Tested Applications</b> | FC                                                                            |
| <b>Preservatives</b>       | PBS with 0.05% Proclin300, 1% rAlbumin                                        |
| <b>Concentration</b>       | 0.5 mg/mL                                                                     |
| <b>Storage</b>             | Store at 2~8°C and protected from prolonged exposure to light. Do not freeze. |
| <b>Note</b>                | For research use only                                                         |
| <b>Isotype</b>             | IgM, $\kappa$                                                                 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                               |