

Product Datasheet

Mouse IgM, κ Isotype Control (FITC) (orb1187573)

Description	Mouse IgM, κ Isotype Control (FITC)
Species/Host	Mouse
Conjugation	FITC
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Concentration	0.5 mg/mL
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	IgM, κ
Expiration Date	12 months from date of receipt.