

Product Datasheet

Hepatitis B Virus Surface Antigen antibody (Biotin) (orb1187600)

Description	Mouse polyclonal antibody against Hepatitis B Virus Surface Antigen conjugated to Biotin
Species/Host	Mouse
Reactivity	Virus
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Hepatitis B Virus Surface Antigen
Preservatives	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide.
Concentration	1mg/mL
Storage	Store at -20°C. Avoid repeated freeze / thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	ELISA: Use at an assay dependent dilution.
Expiration Date	12 months from date of receipt.