

Product Datasheet

GSK3 beta antibody (orb1187636)

Description	Rabbit polyclonal antibody to GSK3 beta
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant protein corresponding to mouse GSK3 beta.
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/mL
Storage	Store at -20°C. Avoid repeated freeze / thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB (1/2000 - 1/5000), IH (1/100 - 1/200)
Expiration Date	12 months from date of receipt.