

Product Datasheet

AP2 alpha/beta antibody (orb1187675)

Description	Rabbit polyclonal antibody to AP2 alpha/beta
Species/Host	Rabbit
Reactivity	Gallus, Human, Mouse, Porcine, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	ChIP, EMSA, IF, IHC, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human AP2 alpha/beta. The exact sequence is proprietary.
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/mL
Storage	Store at -20°C. Avoid repeated freeze / thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P01137
Entrez	7040
Dilution Range	WB (1/500 - 1/1000), IH (1/100 - 1/200), IF/IC (1/100 - 1/500), IP (1/10 - 1/100)
Expiration Date	12 months from date of receipt.