

Product Datasheet

His-tag antibody (Biotin) (orb1187733)

Description	Mouse monoclonal antibody against His-tag conjugated to Biotin
Species/Host	Mouse
Conjugation	Biotin
Tested Applications	ELISA, WB
Immunogen	KLH-conjugated synthetic peptide HHHHHH.
Preservatives	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide.
Concentration	1mg/mL
Storage	Store at -20°C. Avoid repeated freeze / thaw cycles.
Note	For research use only
Clonality	Monoclonal
Dilution Range	ELISA: Use at an assay dependent dilution.
Expiration Date	12 months from date of receipt.