

Product Datasheet

Alpha-actin-1 antibody (orb1187776)

Description	Mouse monoclonal antibody to Alpha-actin-1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence of human Alpha-actin-1. The exact sequence is proprietary.
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/mL
Storage	Store at -20°C. Avoid repeated freeze / thaw cycles.
Note	For research use only
Clonality	Monoclonal
Uniprot ID	P68133
Entrez	58
Dilution Range	WB (1/5000 - 1/10000), IH (1/200 - 1/500), IP (1/100 - 1/200)
Expiration Date	12 months from date of receipt.