

Product Datasheet

Myc-tag antibody (HRP) (orb1187802)

Description	Rabbit polyclonal antibody against Myc-tag conjugated to HRP
Species/Host	Rabbit
Conjugation	HRP
Tested Applications	WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence of Myc-tag-HRP labeled. The exact sequence is proprietary.
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/mL
Storage	Store at -20°C. Avoid repeated freeze / thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB (1/1000 - 1/3000)
Expiration Date	12 months from date of receipt.