

Product Datasheet

GATA1 (pS142) antibody (orb1187932)

Description	Rabbit polyclonal antibody to GATA1 (pS142)
Species/Host	Rabbit
Reactivity	Canine, Human, Monkey, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human GATA1 (pS142). The exact sequence is proprietary.
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/mL
Storage	Store at -20°C. Avoid repeated freeze / thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P15976
Entrez	2623
Dilution Range	WB (1/500 - 1/1000), IH (1/100 - 1/200), IP (1/10 - 1/100)
Expiration Date	12 months from date of receipt.