

Product Datasheet

c-Myc (pT58) antibody (orb1188142)

Description	Rabbit polyclonal antibody to c-Myc (pT58)
Species/Host	Rabbit
Reactivity	Bovine, Canine, Gallus, Human, Mouse, Porcine, Rat, Sheep, Zebrafish
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the N-term region of human c-Myc. The exact sequence is proprietary.
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/mL
Storage	Store at -20°C. Avoid repeated freeze / thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P01106
Entrez	4609
Dilution Range	WB (1/500 - 1/1000), IH (1/100 - 1/200), IP (1/10 - 1/100)
Expiration Date	12 months from date of receipt.