

Product Datasheet

HSC70 antibody (orb1188167)

Description	Rabbit polyclonal antibody to HSC70
Species/Host	Rabbit
Reactivity	Bovine, Gallus, Human, Monkey, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	IF, IP, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human HSC70. The exact sequence is proprietary.
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/mL
Storage	Store at -20°C. Avoid repeated freeze / thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P11142
Entrez	3312
Dilution Range	WB (1/500 - 1/1000), IF/IC (1/100 - 1/500), IP (1/10 - 1/100)
Expiration Date	12 months from date of receipt.