

Product Datasheet

Alpha-synuclein (pY125) antibody (orb1188280)

Description	Rabbit polyclonal antibody to Alpha-synuclein (pY125)
Species/Host	Rabbit
Reactivity	Bovine, Human, Monkey, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human Alpha-synuclein. The exact sequence is proprietary.
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/mL
Storage	Store at -20°C. Avoid repeated freeze / thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P37840
Entrez	6622
Dilution Range	WB (1/500 - 1/2000), IH (1/50 - 1/200)
Expiration Date	12 months from date of receipt.