

Product Datasheet

GATA1 antibody (orb1188435)

Description	Rabbit polyclonal antibody to GATA1
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse, Porcine, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	ChIP, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human GATA1. The exact sequence is proprietary.
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/mL
Storage	Store at -20°C. Avoid repeated freeze / thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P15976
Entrez	2623
Dilution Range	WB (1/500 - 1/1000), ChIP (Use at an assay dependent concentration)
Expiration Date	12 months from date of receipt.