

Product Datasheet

TAU (pS579) antibody (orb1189092)

Description	Rabbit polyclonal antibody to TAU (pS579)
Species/Host	Rabbit
Reactivity	Bovine, Gallus, Human, Mouse, Rabbit, Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human TAU. The exact sequence is proprietary.
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/mL
Storage	Store at -20°C. Avoid repeated freeze / thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P10636
Entrez	4137
Dilution Range	WB (1/500 - 1/1000)
Expiration Date	12 months from date of receipt.