

Product Datasheet

IGF1 Receptor antibody (orb1189110)

Description	Rabbit polyclonal antibody to IGF1 Receptor
Species/Host	Rabbit
Reactivity	Bovine, Gallus, Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human IGF1 Receptor. The exact sequence is proprietary.
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/mL
Storage	Store at -20°C. Avoid repeated freeze / thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P08069
Entrez	3480
Dilution Range	WB (1/500 - 1/1000)
Expiration Date	12 months from date of receipt.