

Product Datasheet

GLUR4 (pS862) antibody (orb1189118)

Description	Rabbit polyclonal antibody to GLUR4 (pS862)
Species/Host	Rabbit
Reactivity	Gallus, Human, Monkey, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human GLUR4. The exact sequence is proprietary.
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/mL
Storage	Store at -20°C. Avoid repeated freeze / thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P48058
Entrez	2893
Dilution Range	WB (1/500 - 1/1000)
Expiration Date	12 months from date of receipt.