

Product Datasheet

JNK3 antibody (orb1189847)

Description	Rabbit polyclonal antibody to JNK3
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IP, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human JNK3. The exact sequence is proprietary.
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/mL
Storage	Store at -20°C. Avoid repeated freeze / thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P53779
Entrez	5602
Dilution Range	WB (1/500 - 1/1000), IF/IC (1/100 - 1/500), IP (1/10 - 1/100)
Expiration Date	12 months from date of receipt.