

Product Datasheet

Estrogen Receptor beta antibody (orb1189921)

Description	Rabbit polyclonal antibody to Estrogen Receptor beta
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant full length protein of human Estrogen Receptor beta. The exact sequence is proprietary.
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/mL
Storage	Store at -20°C. Avoid repeated freeze / thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q92731
Entrez	2100
Dilution Range	WB (1/500 - 1/2000), IH (1/50 - 1/200)
Expiration Date	12 months from date of receipt.