

Product Datasheet

Histone H3 (AcK27) antibody (orb1190001)

Description	Rabbit polyclonal antibody to Histone H3 (AcK27)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A synthetic acetylated peptide corresponding to residues surrounding K27 of human Histone H3 (AcK27)
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/mL
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q16695
Entrez	8290
Dilution Range	WB (1/500 - 1/2000), IP (1/50 - 1/100)
Expiration Date	12 months from date of receipt.