

Product Datasheet

Histone H4 (AcK16) antibody (orb1190625)

Description	Rabbit polyclonal antibody to Histone H4 (AcK16)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence of human Histone H4. The exact sequence is proprietary.
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/mL
Storage	Store at -20°C. Avoid repeated freeze / thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P62805
Entrez	121504
Dilution Range	WB (1/1000 - 1/2000), IH (1/200 - 1/500)
Expiration Date	12 months from date of receipt.