

Product Datasheet

NMDAR2B (pY1336) antibody (orb1191034)

Description	Rabbit polyclonal antibody to NMDAR2B (pY1336)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human NMDAR2B. The exact sequence is proprietary.
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/mL
Storage	Store at -20°C. Avoid repeated freeze / thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q13224
Entrez	2904
Dilution Range	WB (1/500 - 1/1000), IH (1/100 - 1/200)
Expiration Date	12 months from date of receipt.