

Product Datasheet

NDRG3 antibody (orb1193408)

Description	Rabbit polyclonal antibody to NDRG3
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human NDRG3. The exact sequence is proprietary.
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/mL
Storage	Store at -20°C. Avoid repeated freeze / thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9UGV2
Entrez	57446
Dilution Range	WB (1/500 - 1/1000), IH (1/50 - 1/100)
Expiration Date	12 months from date of receipt.